DHERMA ROLLER

RODILLO DE MICROAGUJAS DE TITANIO

TITANIUM MICRONEEDLE ROLLER

NEURONE[®]
COSMÉTICA

DESCRIPCIÓN:

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DESCRIPTION:

Auxiliary to treat hair thinning, contributing to the decrease of hair loss.

- The roller helps remove dead skin from the scalp, helping to stimulate the skin to create the right environment for the application of Capissen Lotion.
- Contributes to optimal and healthy hair growth.
- Contributes to collagen production and cell renewal.
- By combining a scalp treatment with Neurone Dherma Roller, it helps for more effective absorption.

DIRECTIONS:

- Before first use clean your Neurone Dherma Roller. The correct way to clean your Roller is to immerse the rotating head in 75% alcohol for 5 to 8 minutes. Never use your roller without first cleaning it.
- On the scalp, apply Capissen Lotion and massage gently.
- Roll the roller in one direction with gentle pressure on the scalp, starting at the front of the head at the hairline and working towards the back.
- Go over the area 5 to 10 times changing direction. The roller must be used in all directions (vertical, horizontal and diagonal).
- Use only once a day, maximum 3 times a week.

TYPE OF SKIN:

- Do not use if you have sensitive skin or are allergic to metals.
- The roller should not be used on open wounds, active acne, pustular, eczema, sunburn or abnormal skin diseases.

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PRECAUTIONS:

- Do not reuse the head without cleaning it before next use.
- Do not share NEURONE DHERMA ROLLER with other people.
- Do not ride more than indicated in the instructions.
- If your skin shows hypersensitive responses when treated, consult your doctor.
- Store the roller in its box after each use.
- Keep out of the reach of children.